



# A capillary furnace designed for variable-temperature X-ray diffraction

**Branden Lowe<sup>1</sup>, William O'Donnell<sup>2</sup>, Selva Venilla Raju<sup>2</sup> and Ravhi Kumar<sup>2</sup>**

<sup>1</sup> Whitworth University, Spokane, Washington 99251

<sup>2</sup> University of Nevada Las Vegas, Las Vegas, Nevada 89154




## Introduction

In the present work, we build an inexpensive capillary furnace to be used in an X-ray diffraction machine. The basis of the furnace comes from work done by Lavigueur and colleagues. The purpose of this project is two-fold, namely, design a furnace that is inexpensive and conduct variable-temperature X-ray diffraction on air-stable samples. Capillary temperatures can reach temperatures up to 1000 K. We attach the furnace to a standard goniometer head, which is then mounted in the X-ray diffraction machine. The simplicity inherent to the furnace allows for easy building and operation.

## Objectives

▪ **Design a furnace that is capable of reaching 1000 K and can be installed in an X-ray diffraction machine.**

▪ **X-ray diffraction machine:**



▪ **Inexpensive: under \$1,000**

▪ **Furnace quote from a company: \$68,800**

▪ **More advanced furnace**

▪ **Easy operation:**

▪ **Loading capillaries**

▪ **Assembling furnace**

▪ **Replacing worn heating elements**

▪ **Vacuum capability**

▪ **Next step in furnace advancement**

## Methods and Materials

▪ **Solidworks 2011 was used to design the furnace**

▪ **MACOR ceramic used to make furnace**

▪ **Max temperature: 1000° C**



▪ **Sustainable at 800° C**

▪ **Heating Element: Nickel Chromium (Nichrome) wire**

▪ **Melting Temperature: 1350° C**

▪ **Two models were created:**

▪ **Model A and Model B:**



▪ **Why we chose model B:**

▪ **Reduce chances of short circuiting**

▪ **Improved quartz casing seal – minimize heat loss**

## Prototype

▪ **Goals:**

▪ **Verifying the literature**

▪ **Testing various gauges of Nichrome wire**

▪ **Tested: 20, 24, 26, and 28 gauge wire – chose 24**

▪ **24 gauge offered best temperature control, wire controllability, and amperage usage.**

▪ **Design:**




▪ **Quartz tube used to withstand temperatures**

- **Ceramic insulators to reduce chances of wood igniting**
- **Shortcomings:**
  - **Crimping device finicky – resistance varied**
  - **Poor seal between quartz and wood – loss of heat after thermal cycling**
  - **Melting temperature of Nichrome**
    - **In the future use platinum or tungsten**


## Final Design

▪ **Furnace unassembled:**



To left, 5 ceramic pieces, goniometer head, quartz tube, 2 threaded bolts and nichrome wire

▪ **Furnace Assembled:**



## Acknowledgements

I would personally like to thank Paul Ellison for machining the heater, Daniel Antonio and Jason Baker for familiarizing us with both the equipment and the facilities here, Professors Andrew Cornelius, John Farley, and UNLV for giving me the opportunity, and lastly, the REU program of the National Science Foundation under grant DMR-1005247 is gratefully acknowledged.

## References

Christine Lavigueur, E. Johan Foster and Vance E. Williams. A simple and inexpensive capillary furnace for variable- temperature X-ray Diffraction. *Applied Crystallography*, November 2007

**Funding Source: The National Science Foundation**